

Lecture 4 – Vectors for Constructing Genomic and cDNA Libraries

Cloning Vectors

- **Purpose**

- ✓ **Facilitates** the insertion, storage and manipulation of exogenous DNA

- **Types**

- ✓ **Plasmids**
- ✓ **Bacteriophages**
- ✓ **Cosmids**
- ✓ **BACs, PACs and YACs**


Plasmids

- Plasmid vectors are **double-stranded**, circular, self-replicating, **extra-chromosomal** DNA molecules.
- **Advantages:**
 - ✓ Small, easy to handle
 - ✓ Straightforward selection strategies
 - ✓ Useful for cloning small DNA fragments (< 10kbp)
 - ✓ **Replicate autonomously**
- **Disadvantages:**
 - ✓ Less useful for cloning large DNA fragments (> 10kbp)


Features of Plasmid Vectors

- Contain an **origin of replication**, allowing for replication independent of host's genome
- Contain **selective markers** for antibiotic resistance and blue-white screening
- Contain a **multiple cloning site (MCS)**
- **Easy** to isolate from host cell

Examples of Plasmid Vectors



- pUC19 is one of a series of plasmid cloning vectors created by Joachim Messing et al. at UC Davis.
- The designation "pUC" is derived from the classical "p" prefix (denoting "plasmid") and the abbreviation for the University of California, where early work on the plasmid series had been conducted.



- pBR322 was one of the first widely used *E. coli* cloning vectors created in Herbert Boyer laboratory at UC San Francisco.
- It was named after the postdoctoral researchers who constructed it: Bolivar and Rodriguez.

Bacteriophages

Lambda (λ) phage

- Lambda (λ) phage is a virus particle consisting of a head, containing **double-stranded linear DNA (ds DNA)** as its genetic material, and a tail that can have tail fibers.
- The phage particle recognizes and binds to its host, *E. coli*, causing DNA in the head of the phage to be **ejected** through the tail into the cytoplasm of the bacterial cell.
- Usually, a “**lytic cycle**” ensues, where the lambda DNA is replicated many times and the genes for head, tail and lysis proteins are expressed.
- This leads to assembly of **multiple new** phage particles within the cell and subsequent cell lysis, releasing the cell contents, including **virions** that have been assembled, into the environment.

M13 phage

- M13 is a **filamentous** bacteriophage composed of **circular single stranded DNA (ssDNA)**, which is 6407 nucleotides long encapsulated in approximately 2700 copies of the major coat protein P8, and capped with 5 copies of two different minor coat proteins (P9, P6, P3) on the ends.

Bacteriophages



- **Advantages:**

- ✓ Useful for cloning large DNA fragments (10-23 kbp)
- ✓ Inherent size selection for large inserts


- **Disadvantages:**

- ✓ Less easy to handle than plasmid vectors


λ Phage




Life Cycle of λ Phage



Lysis Plaques of Lambda Phage on *E. coli*



λ Phage Vector



Lambda Vector

- Infectious λ can be assembled *in vitro*.
- Foreign DNA can be incorporated into the λ genome.
- Non-essential genes are removed.
- Phage assembly can occur with 40-52 kb of DNA



- accommodates up to 7-10 kb of foreign DNA (depending on the vector)
- 13 kb “stuffer fragment” (lysogeny genes) discarded
- accommodates 11-20 kb of foreign DNA


Lambda Insertion Vectors

- Insertion vectors are the **simplest** form of lambda cloning vectors.
- The vector itself can be grown (must contain at least 75% of the wild-type genome length).
- Foreign DNA fragments are inserted into **unique** restriction sites in the vector genome.
- Packaging requirements **limit** insert fragment size to 0-10 kb due to the limitations on viral genome size (75% to 105% of the wild-type length = 50 kb).
- In the construction of insertion vectors, 25% of the dispensable portion of the wild-type lambda genome is **eliminated**.
- Because the lambda genome approximates a random 50-kb DNA sequence, most 6-mer restriction enzymes will cut **several times** (average frequency is 1/4000 bp).
- Removal of **unwanted** restriction sites involves **mutagenesis** of the phage DNA followed by **selection** for efficient growth on hosts carrying **different** restriction/modification systems.
- After **multiple** rounds of mutagenesis and selection, phage genomes lacking unwanted restriction sites are **recovered**.


Lambda Replacement Vectors

- Lambda replacement vectors have a loading capacity of 10-23 kb and consists of a **full-length** lambda molecule with two **identical** restriction sites flanking a non-essential region called a **stuffer** fragment.
- A stuffer fragment is deleted and replaced by foreign DNA.
- A **selection system** is required to differentiate between wild-type and recombinant phage.
- This is done by placing **relevant** genes onto the stuffer fragment, the **loss** of which gives rise to a **detectable** phenotypic signal.
- One example is:
 - ✓ LacZ inactivation by inserting the lacZ gene onto the stuffer fragment
 - ✓ Loss of the fragment and, therefore, lacZ can be monitored by the well known blue/white color discrimination.


Constructing a Genomic Library in λ Phage



M13 Phage




Lysis Plaques of M13 Phage on *E. coli*



- M13 and the other filamentous DNA phages are unusual in that they are released from the infected cell, **without** causing the death of the bacterial cell.
- Therefore, phage plaques of M13 are seen as **areas of reduced growth** in the bacterial lawn, i.e., clear zones of bacterial lysis do not occur.

M13 Vector



Bacteriophage M13

- Bacteriophage M13 has a single-stranded DNA genome.
- After M13 infection into a bacterial cell, a complementary DNA strand is synthesized, generating a double-stranded **replicative form (RF)** of the bacteriophage genome.
- The **complementary** strand then serves as a **template** for synthesis of **new** single-stranded viral DNA by a rolling circle mechanism, except that the second strand is not synthesized until a new cycle of replication is initiated.
- The single-stranded DNA is cut into **genome-length** fragments and extruded from the cell.
- This single-stranded DNA is useful in sequencing studies on foreign DNA cloned into M13.
- In addition, M13 clones can easily be subjected to **site-directed** mutagenesis.

Cloning Into an M13 Vector

- M13 vectors have been engineered to contain restriction endonuclease cut sites in the double-stranded **replicative** form.
- For cloning, the replicative form can be isolated from infected bacteria or generated **artificially** from the single-stranded form with a complementary synthetic oligonucleotide primer and DNA polymerase plus ligase.
- The double-stranded form is cut with a restriction endonuclease and the DNA is inserted much like any other cloning procedure.
- The double-stranded replicative form with the insert is infected into bacteria and generates **single-stranded clones** by the normal process of single-stranded genomic replication.

Cosmids

- **Features**


- ✓ A cosmid is a type of hybrid plasmid that contains a lambda phage **cos sequence**, which is the sequence that is cut to produce the cohesive, single-stranded extensions located at the ends of the linear DNA molecules of certain phages such as lambda.
- ✓ Cosmid DNA sequences originate from the lambda phage (**cos sites + plasmid = cosmid**).
- ✓ Cos sequences are **~200 base pairs** long and **essential** for packaging.
- ✓ Cosmids are **recombinant plasmids** contained in phage shells.

- **Advantages**

- ✓ Useful for cloning very large DNA fragments
(32 - 47 kbp)
- ✓ Inherent size selection for large inserts
- ✓ Handle like plasmids

Cosmid Engineering


- Cosmids combine essential elements of a plasmid and lambda systems.
- Fragments from 30 - 46 kb can be accommodated by a cosmid vector.
- Cosmids are extracted from bacteria and mixed with restriction endonucleases.
- Cleaved cosmids are mixed with foreign DNA that has been cleaved with the same endonuclease.
- Recombinant cosmids are packaged into lambda capsids.
- Recombinant cosmid is injected into the bacterial cell where the R cosmid arranges into a circle and replicates as a plasmid.
- It can be maintained and recovered just as plasmids.




- Shown above is a 50,000 base-pair long DNA molecule bound with six EcoRI molecules, and imaged using an atomic force microscope.
- This image clearly indicates the six EcoRI "sites" and allows an accurate restriction enzyme map of the cosmid to be generated.

<http://homer.ornl.gov/cbps/afmi/maging.htm>

Example of Cosmid Vector



Constructing a Genomic Library in Cosmids



BACs, PACs and YACs

- Features

- ✓ BAC: Bacterial Artificial Chromosome – DNA construct based on a **functional fertility plasmid (F-plasmid)**
- ✓ PAC: P1-derived Artificial Chromosome – DNA construct derived from DNA of **P1 bacteriophage**
- ✓ YAC: Yeast Artificial Chromosome – originated from a **bacterial plasmid**; contains additionally a yeast **centromeric region (CEN)**; **autonomously replicating sequence (ARS)** which is a yeast origin of DNA replication; a cluster of **unique** restriction sites; and a **selectable marker** and a **telomere region** at the end of each arm


- Advantages

- ✓ Useful for cloning **extremely large** DNA fragments (100 - 2,000 kbp), depending on the artificial chromosome vector used

- Disadvantages


- ✓ **Not easy** to handle extremely large DNA molecules

Example of a BAC




- Cre-Lox recombination consists of a single enzyme, Cre recombinase, which recombines a pair of short target sequences called the *Lox* sequences.
- This system can be implemented without inserting any extra supporting proteins or sequences.
- The Cre enzyme and the original *Lox* site called the *LoxP* sequence are derived from bacteriophage P1.


Cloning in a BAC



Example of a PAC




Example of a YAC




***Capable of carrying inserts of 200-2000 kbp in yeast**

Cloning in a YAC



Constructing a Genomic Library in a YAC



ARS – autonomously replicating sequence
CEN – yeast centromeric region

Expression Vectors

- Expression vectors enable a particular gene to be not only **amplified** but also **expressed** in a host cell.
- The vector is constructed to contain appropriate **regulatory** sequences, such as a promoter and operator, so that the host-cell machinery can transcribe the gene and translate the resultant messenger RNA to synthesize the corresponding protein.
- Such vectors are thus **essential** for the manufacture of, for example, mammalian proteins by bacterial host cells.
- Expression vectors used in prokaryotes are typically based on **plasmids** or **phages**, or **plasmid–phage hybrids** (**phagemids**).
- Some **eukaryote** proteins are **extensively modified** during or following their synthesis, for example by the addition of carbohydrate groups.
- **Prokaryote** host cells are **unable** to accomplish these modifications, and therefore expression systems based on eukaryotic cells must be used instead.
- For example, vectors based on baculovirus, a DNA virus of insects, are used successfully in cultures of insect cells.
- **Secretion vectors** allow both **expression** and **secretion** of the novel protein by the host cell, by ensuring that the expressed protein carries a **signal peptide** that allows it to be **transported** across the plasma membrane.