

Lecture 2: How to Clone a Gene

Cloning

- Cloning utilizes **specialized** DNA technology to produce **multiple, exact copies** of a **single** gene or other segment of DNA to obtain enough material for further study.
- The resulting cloned (copied) collections of DNA molecules are called **clone libraries**.
- A second type of cloning exploits the **natural process** of cell division to make many copies of an **entire** cell.
- The genetic makeup of these cloned cells, called a **cell line**, is identical to the original cell.
- A third type of cloning produces **complete, genetically identical** animals, such as the famous Scottish sheep, Dolly, and plants such as rice.
- Often times, making **subclones** will be necessary, which involves breaking a cloned fragment into **smaller** fragments for **further** cloning.

Cloning a gene involves making a recombinant DNA molecule.

- A recombinant DNA molecule is one that **combines** DNA from two sources to create **new** genetic material, which is able to **replicate** and **express** under certain conditions.
- **Stanley Cohen** and **Herbert Boyer** created the **first** recombinant DNA **organism** (1973) using recombinant DNA techniques pioneered several years earlier by **Paul Berg**.
- An example of a recombinant is the **human gene for insulin** which was cloned into a carrier DNA (vector).
- The resulting circular DNA molecule was placed in **bacteria** to produce insulin in **large** quantities for treating diabetics.
- Such bacteria are recombinant organisms.
- It is possible to create genetically modified organisms because of the **universal nature** of the genetic code.

Basic Components of Recombinant DNA Technology

- Restriction endonucleases – enzymes that **cleave** DNA into pieces
- Insert – **target** DNA fragment or gene
- Vector – an **engineered** DNA molecule used to **transfer** and **propagate** various insert DNAs
- DNA ligase – enzyme that **joins** two separate DNA molecules together at the **ends**
- Host – **surrogate** cell harboring recombinant DNA for **amplification** and possible **expression**
- Techniques to **deliver** recombinant DNA into the host

Nucleases

- **Endonuclease** – nucleases that cleave phosphoester bonds **within** a nucleic acid chain
 - They may be **specific** for **RNA** or for **single-stranded** or **double-stranded** DNA.
- **Exonuclease** – nucleases that cleave phosphoester bonds **one at a time** from the **end** of a polynucleotide chain
 - They may be **specific** for either the **5' or 3' end** of DNA or RNA.

Restriction Endonucleases


- **Werner Arber, Hamilton Smith and Daniel Nathans** discovered that enzymes (endonucleases) isolated from bacteria cut foreign DNA at precise locations.
- There are **three** types of restriction endonucleases: I, II, and III.
- Type I cleaves **1000 bp away** from the recognition site.
- Type II cleaves DNA at **defined** positions **close to or within** their recognition sequences and produces **discrete** restriction fragments and **distinct** gel banding patterns. Type II enzymes are used in the laboratory for routine DNA analysis and gene cloning.
- Type III cleaves **24-26 bp away** from the site.

Type II Restriction Endonucleases


- More than **3000** type II restriction endonucleases have been discovered.
- They recognize short, usually **palindromic**, sequences of **4–8 bp** and, in the presence of **Mg²⁺**, cleave the DNA within or in close proximity to the **recognition sequence**.
- The precise mechanism of cleavage has **not** yet been established for any restriction enzyme—the main uncertainty concerns the **number** of Mg²⁺ ions directly involved in cleavage.
- Cleavage in the two strands usually occurs in a **concerted fashion** and leads to **inversion** of configuration at the **phosphorus** group.
- The **products** of the reaction are DNA fragments with a **3'-OH** and a **5'-phosphate**.

Restriction Endonucleases

Staggered cut




Blunt cut




Restriction endonucleases can be used to cleave DNA into **defined** fragments.

Creating Recombinant DNA Molecules




Gel Electrophoresis of DNA Cut With Different Restriction Enzymes



Restriction Map of Recombinant Plasmid Vector Created by GenScript

- GenScript pDream2.1/MCS is an excellent expression vector.
- There are **seven** restriction enzyme sites in MCS.
- A gene cloned into MCS can be expressed in any one of the **three** major protein expression systems: bacteria, Insect cells and mammalian cells.



Flag Tag → *BamH I*

1804 TAT ACC **ATG** GAT TAC AAG GAT GAC GAC GAT AAG GGA TCC
 ATA TGG TAC CTAATG TTC CTA CTG CTG CTA TTC CCT AGG

Met Asp Tyr Lys Asp Asp Asp Asp Lys Gly Ser

Bgl II *EcoR I* *Kpn I* *Sac I* *Hind III* *Not I*

1843 GAG ATC TG **AAT TCG GTA CCG AGC TCA AGC TTG CGG CCG**
 CTC TAG ACC **TTAAGC CAT GGC TCG AGT TCG AAC GCC GGC**


Glu Ile Trp Asn Ser Val Pro Ser Ser Ser Leu Arg Pro

SP6 Primer

1882 CCTATAGTGT**CACCTA**AATCGTAACCCAGC
 GGATATCACAGTGGATTTAGCATTGGGTCG

Restriction Map




A map can be generated by using the **overlaps** between the fragments generated by **different** restriction enzymes.



A restriction map is a **linear sequence** of sites separated by **defined** distances on DNA.

Vectors


- Vectors are DNA molecules that can be **moved into** and **replicated** in an organism **independent** of the host.
- They are shuttles for **moving** and **copying** another DNA fragment or gene.
- Vectors are classified by:
 - ✓ the **organisms** which replicate the vector
 - ✓ the **size** of DNA that can be inserted
- A vector contains an **origin of replication (ori)**, **selective marker(s)** and **multiple cloning site (MCS)**.



Common Vectors

Vector	Host	Size of insert accepted (kb)
Plasmid	<i>E. coli</i>	<15
Cosmid (plasmid/ λ site)	<i>E. coli</i>	<50
Fosmid (F plasmid ORI + λ cos site)	<i>E. coli</i>	<50
Bacteriophage	<i>E. coli</i>	<100
P1(contains deletion in λ)	<i>E. coli</i>	<125
Bacterial artificial chromosome (BAC)	<i>E. coli</i>	100-500
Yeast artificial chromosome (YAC)	Yeast	250-1,000
P1-derived Artificial Chromosome (PAC)	<i>E. coli</i>	70-100 varies
Ti plasmid	Plant	
Viral vector (retrovirus, lentivirus, adenovirus, adeno-associated virus)	Mammalian cells	varies varies
Engineered nanoparticles	Mammalian cells	

Plasmid – Circular Extrachromosomal DNA





Cloning vectors can be specialized for different purposes.

- Cloning a fragment of DNA requires **specially engineered** vectors.
- **Shuttle** vectors can be propagated in **more** than one type of host cell.
- **Expression** vectors contain **promoters** that allow transcription of any cloned gene.
- They generally contain a **multiple cloning site** (MCS)—a sequence of DNA containing a series of tandem restriction endonuclease sites.


DNA Ligase

DNA ligase is an enzyme, identified by **Martin Gellert** in 1967, that **joins** two DNA molecules by forming two covalent **phosphodiester bonds** between 3' hydroxyl ends of one nucleotide ("acceptor") with the 5' phosphate end of another ("donor").




Cloning: plasmid transformation

- A plasmid contains **three** key sites: **ori**, **antibiotic-resistant gene** and **lacZ** with an **MCS**, together with the insert DNA to be cloned and an appropriate restriction enzyme.
- Restricted insert fragment and vector will be **combined** and **ligated** together.
- The final pool of this DNA can be **transformed** into *E. coli*.



Transformation – acquisition of new genetic material by incorporation of added exogenous, nonviral DNA



Selecting Recombinant Molecules


Three types of cells can result from attempts to introduce a DNA molecule into a bacterial cell:

- Cells **without** vector
- Cells **with** vector but **no** inserted DNA
- Cells **with** vector and **inserted** DNA

Selecting for Cells with Vector

- Vectors are commonly engineered to carry **antibiotic-resistance** genes.
- Host bacteria **die** in the presence of the antibiotic.
- Bacteria harboring the vector **survive**.
- Growing cells on media with antibiotics **ensures** that all growing cells **must** carry the vector.


Selecting for Cells with Vector




Selecting for Cells with Inserted DNA

- The **site of insertion** of the DNA of interest can be within a gene on the vector.
- Insertion of a DNA fragment will **disrupt** the vector gene.
- The vector gene *lacZ* produces an enzyme which allows the bacteria to turn **blue** in the presence of **certain** media.
- Insertions in the *lacZ* gene **prevent** *lacZ* enzyme production and the bacteria are **white**.
- Bacteria with vector that are white carry a DNA **inserted** in the *lacZ* gene. (Rare mutations in the *lacZ* gene also will render white bacteria.)

Selecting for Cells with Inserted DNA




Recombinant Selection with pUC8



Blue-white Colony Selection

Blue/white selection allows the **identification** of bacteria that contain the **vector plasmid** and vector plasmids that contain an **insert**.

E. coli colonies on agar plates with ampicillin, IPTG, and the color indicator X-gal



Blue/white Colony Screening Assay, Distinguishing Recombinant Colonies (white) Among Non-recombinant Ones (blue)

- X-Gal (5-bromo-4-chloro-3-indolyl-beta-D-galactopyranoside) is an inert chromogenic substrate for beta-galactosidase, an enzyme that promotes lactose utilization.
- Beta-galactosidase hydrolyzes X-Gal into a colorless galactose and 4-chloro-3-brom-indigo which forms an intense blue precipitate. Induction of the *lacZ* gene with IPTG leads to the hydrolysis of X-Gal and to the development of blue colonies.
- The IPTG is a highly stable synthetic analog of lactose. It inactivates the *lac* repressor and induces synthesis of beta-galactosidase, an enzyme that promotes lactose utilization.



Delivery of Recombinant DNA Molecules


Introduction methods focus on getting recombinant DNA molecule (vector and insert) across the membrane of the host cell:

- ✓ **Electroporation** – electricity opens a temporary hole
- ✓ **Microinjection** – microscopic needle injects DNA
- ✓ **Liposomes** – fatty bubbles move across membranes
- ✓ **Particle bombardment** – DNA coated “bullets” are shot into cell
- ✓ **Chemicals** – salts open a temporary hole
- ✓ **Viruses** – DNA or RNA viruses infect cell
- ✓ ***Ti* plasmid** – tumor-inducing genes are removed, gene of interest inserted and modified plasmid introduced into plant protoplasts

Applications of Recombinant DNA Technology

- **Large-scale production of proteins by genetically engineered bacteria.**
 - ✓ **Examples: insulin, growth hormone, interferon's and blood clotting factors (VIII and IX).**
- **Vaccine development**
- **Gene therapy for genetic diseases**
- **Food production**
- **Stem cell research: iPS (induced pluripotent stem cells)**

Recombinant DNA Technology



Commercial Fermenters

