

Lecture 8: Analysis of Gene Expression (Part 1, Global Gene Expression)

Global Gene Expression


- A **new** field of genetic analysis called global gene expression has emerged in recent years, driven by the realization that traditional techniques of **linkage** and **association analysis** can be applied to **thousands** of transcript levels measured by **microarrays**.
- Genetic dissection of **transcript abundance** has shed light on the **architecture of quantitative traits**, provided a new approach for connecting **DNA sequence variation with phenotypic variation** and improved our understanding of **transcriptional regulation** and **regulatory variation**.

Global Gene Expression

- Studies in a number of species have documented abundant **heritable variation** in gene expression among individuals and strains.
- Linkage analysis has been used to map **thousands of loci** that affect gene expression.
- Expression traits consistently show **complex inheritance**, explicable only by multiple underlying loci and possibly interactions among the loci.
- Many types of **genetic complexity** are observed across the thousands of **expression traits**.
- A locus that **affects** gene expression can be classified according to its **location** as “**local**” (**near** the genomic location of the gene) or “**distant**” (**elsewhere** in the genome).
- Many gene expression traits are affected by **local** regulatory variation.
- It appears that **most**, but not all, local regulatory variations function in **cis**, with perhaps a **quarter to a third** acting in **trans**.
- Most transcripts link to loci that are **distant** from the genomic locations of the genes that encode the **corresponding** transcripts.
- One **common** feature of these loci is “**hot spots**”: **individual** loci that affect **large** numbers of transcripts.
- Future research will focus on **high-throughput** identification and characterization of polymorphisms that affect expression, extending studies to population samples and applying the approach to other **global molecular phenotypes**.

Heritability (H)

H estimates the proportion of the phenotypic variation in a population due to genetic differences.



- Genetic variants are determined mostly by the **additive** effects of **recessive** alleles of different genes.
- Variants also can be influenced by the effects of a few **dominant** alleles and by **epistasis**.

Cis-acting Elements and Trans-acting Factors

- The **transcription initiation complex** is composed of **promoter sequences** and **DNA binding proteins**.
- These two components of transcription are normally described as **cis-acting elements** and **trans-acting factors**.
- Cis-acting **elements** are DNA sequences in the **vicinity** of the **structural** portion of a gene that are **required** for gene expression.
- Trans-acting **factors** usually are considered to be proteins that **bind** to the cis-acting sequences to **control** gene expression.

Human Transcription Complex

Activators


These regulatory proteins bind to DNA at distant sites known as enhancers. When DNA folds so that the enhancer is brought into proximity with the transcription complex, the activator proteins interact with the complex to increase the rate of transcription.

Repressors - DNA binding proteins

These regulatory proteins bind to "silencer" sites on the DNA, preventing the binding of activators to nearby enhancers and so slowing transcription.

Basal factors

These transcription factors, in response to coactivators, position RNA polymerase at the start of a protein-coding sequence, and then release the polymerase to transcribe the mRNA.



Coactivators


These transcription factors transmit signals from activator proteins to the basal factors.

Core promoter

Analysis of Transcription

Two general approaches:

1. Measure the **amount** of a **specific** transcript produced by a cell or a group of cells.
2. Measure the transcriptional **activity** of a particular gene or genes under certain conditions.



Cautions for Handling RNA

Preparation for northern analysis:

1. **EVERYTHING** must be thoroughly cleaned of **RNAases**, and gloves must be worn at all times.
2. Use **RNAse-ZAP** or some other commercially available RNase cleaner to clean the work area and all equipment that may come into contact, directly or indirectly, with the RNA.
3. If no cleaners are available, use a dilute solution of **bleach**.
4. Clean everything with **95% ethanol**, which will facilitate **evaporation** of any residual RNase-ZAP.
5. Place **signs** around the work area that stress the importance of **no contamination**.
6. Use **only DEPC-treated** (diethylpyrocarbonate) chemicals.
7. Chemicals can be DEPC-treated by adding two ml of DEPC per liter of solution, followed by **autoclaving**.


Protocol for RNA Isolation by TRIzol and RNeasy Column Binding for Affymetrix Small Sample Preparation

- The basis for this protocol is that **TRIzol** is used to initially **lyse** the cells and **solubilize** nucleic acids and proteins and other biomolecules.
- The organic mixture (**plus chloroform**) is separated using a phase-lock gel tube, which ensures that virtually all of the aqueous phase is collected **without contamination** from the interface or organic layer.
- RNA is then purified from the aqueous phase by **RNeasy** (Qiagen) purification.
- This procedure is preferred over alcohol precipitation for smaller amounts of RNA because of the risk of losing or not seeing a miniscule nucleic acid pellet.
- It also reduces the risk of organic solvent contamination of the RNA solution.
- **Any** of the RNeasy-style technologies can be used : the RNeasy mini-extraction kit, MiniElute kit or the Affymetrix sample clean-up kit.

Quantification of Transcripts

- Northern Blot Analysis
- RNase Protection Assay
- Primer Extension Assay
- Quantitative PCR
- Microarray
- Serial Analysis of Gene Expression
- Differential Display
- Subtractive Hybridization
- RNA sequencing

Northern Blot Analysis




- It is the **gold standard** for quantifying RNA molecules.
- It facilitates **direct** measurement of RNA molecules from samples.
- It provides **information** about a RNA transcript, e.g., size and abundance, alternative splicing variants, etc.)

RNase Protection Assay

- RNase protection assay is used to **identify** individual RNA molecules in a heterogeneous RNA sample extracted from cells.
- The extracted RNA is first mixed with **anti-sense** RNA probes that are complementary to the sequence or sequences of interest and the complementary strands are **hybridized** to form double-stranded RNA.
- The mixture is then exposed to **ribonucleases** that specifically cleave **only** single-stranded RNA but have **no** activity against double-stranded RNA.
- When the reaction runs to completion, susceptible RNA regions are degraded to very short **oligomers** or to individual **nucleotides**.
- The **surviving** RNA fragments are those that are complementary to the added anti-sense strand and contain the **sequence of interest**.
- **Separation** of the hybridization mixture on a polyacrylamide gel and exposure to an X-ray film will **identify** the RNA of interest.
- This technique can identify **one or more** RNA molecules of known sequence at **low** concentration.
- RNase protection assay is **more** sensitive and quantitative than northern blotting and is useful to **map** introns and 5' and 3' ends of transcribed gene regions.



Antisense RNA

- **Antisense** transcripts represent a class of long **noncoding RNAs** (**ncRNAs**).
- Antisense transcripts are transcribed from the strand **opposite** to that of the sense transcript of either protein-coding or non-protein-coding genes.
- RNA can form **duplexes** with its complementary molecule.
- Cells contain genes that are naturally translated into antisense RNA molecules capable of **blocking** the translation of other genes in the cell.
- Apparently, this is another method of **regulating** gene expression.
- Antisense transcripts preferentially accumulate in the nucleus.




cRNA = synthetic transcripts of a specific DNA molecule or fragment, made by an *in vitro* transcription system; RNA derived from cDNA through standard RNA synthesis

RNase Protection Assay



Primer Extension Assay

- In primer extension, the **primer** introduced to the mRNA pool will **hybridize** with the RNA of interest if it is present.
- Hybrids are then extended by reverse transcriptase.
- The information gained through this method includes the confirmation of the **presence** of the RNA of interest, the **location** of the **transcription start site**, and, if an excess of primer is used, the **amount** of the **particular** RNA present in the sample.




Analysis by Microarrays

- A microarray is a tool for analyzing **gene expression** that consists of a small **membrane** or **glass slide** containing **samples** of many genes arranged in a **regular pattern**.
- Microarrays are a significant advance both because they may contain a **very large** number of genes and because of their **small size**.
- Microarrays are useful when one wants to **survey** a large number of genes **quickly** or when the sample to be studied **is small**.
- Microarrays may be used to assay gene expression within a **single sample** or to compare gene expression in **two different** cell types or tissue samples, such as in healthy and diseased tissue.
- Because a microarray can be used to examine the expression of hundreds or thousands of genes at once, it has **revolutionized** the way scientists **examine** gene expression.

Types of Microarrays

- There are **three** basic types of samples that can be used to construct DNA microarrays, **two** are **genomic** and the other is "**transcriptomic**", i.e., it measures **mRNA** levels.
- What makes them **different** from each other is the **kind of immobilized DNA** used to generate the array and, ultimately, the **kind of information** that is derived from the chip.
- The **target DNA** used will also determine the type of control and sample DNA that is used in the hybridization solution.

Spotting a Glass Slide for Microarray Analysis



Colors of a Microarray



Serial Analysis of Gene Expression (SAGE)

- SAGE is used to produce a **snapshot** of the **mRNA population** in a sample of interest in the form of **small tags** that correspond to **fragments** of those transcripts.
- SAGE experiments include the **following steps**:
 - ✓ Isolate mRNA from tissue.
 - ✓ Use a retrotranscriptase to copy mRNA to cDNA.
 - ✓ Extract a small amount of sequence from a defined position of each cDNA molecule.
 - ✓ Link these small pieces of sequence together to form a long chain (or concatemer).
 - ✓ Example of a concatemer:
CATGACCCACGAGCAGGGTACGATGATAC**CATG**GAAACCTATGCACCTTGGGTA
GCAC**CATG**
 - ✓ Clone these chains into a vector which can be taken up by bacteria.
 - ✓ Sequence the chains.
 - ✓ Process this data with a computer to count the small sequence tags.

Features of SAGE


- Provides researchers with a means to study **expression** of an **entire** genome rather than single genes.
- Facilitates **comparison** of gene expression differences between **two** samples with information about **abundance** of transcripts.
- Affords **quantitative** gene expression through library screening.

Two Principles of SAGE

1. **A 9-bp tag** is sufficient to unambiguously identify a gene.
2. **Concatenation** (linking together) of these short DNA sequences increases the efficiency of identifying unique transcripts in a serial manner.

SAGE Protocol

1. mRNA → cDNA
2. Cleave with anchoring enzyme (AE), such as *NlaIII*.
3. Isolate 3'-most transcript of each cDNA by binding to streptavidin beads
4. Divide cDNA in half.
5. Ligate to one of two linkers, each with a tagging enzyme (TE) site.
6. Ligate the two pools of tags together.
7. Ligated linkers serve as primers for amplification.
8. Cleave PCR products with AE to isolate ditags.
9. Concatenate by ligation.
10. Clone.
11. Sequence.





Measurement of Transcriptional Activity

- **Dnase I Footprint Assay**
- **Gel Retardation or Electrophoretic Mobility Shift Assay (EMSA)**
- **Yeast One-hybrid Assay**
- **Reporter Assay**
- **Mutagenesis of DNA Sequences**
- **Chromatin immunoprecipitation (ChIP)**
 - ✓ **ChIP**
 - ✓ **ChIP-chip and ChIP-sequencing**
 - ✓ **Chromosome conformation capture (3C)**

DNase I Footprinting

- DNase I footprinting is a method of studying **DNA-protein** interactions and identifying the DNA sequence to which a protein binds.
- A target DNA fragment about **100-300 bp** in length is either PCR-generated or cut from a vector and then **uniquely** labeled (**at only one end**) and incubated with protein (**usually nuclear extract**), followed by **controlled digestion** with DNase I, which cuts the probe randomly but only **once**.
- The digested DNA is recovered from the reaction and resolved on a polyacrylamide gel along with **G+A chemical sequencing reaction** (Maxam–Gilbert sequencing based on nucleobase-specific partial chemical modification of DNA and subsequent cleavage of the DNA backbone at sites adjacent to the modified nucleotides), which uses the same probe as the template.
- The regions **bound** by proteins will be **protected** from DNase I digestion and will be shown as a **blank** area on the gel track, whereas the **exact** protein-bound sequence can be read out by comparing the **location** of the blank with the sequencing reaction.

DNase I Footprinting




Li et al. *American Journal of Physiology - Gastrointestinal and Liver Physiology* 276:G1425-G1434. 1999.

Gel Retardation Assay

- The gel retardation or **electrophoretic mobility shift** (EMSA, or gel shift) assay is a sensitive technique for studying **protein–DNA** interactions.
- It is widely used to investigate the **formation** of protein-DNA complexes in a **crude** nuclear protein extract.
- The method relies on the **stability** of protein DNA complexes when subjected to **non-denaturing** polyacrylamide gel electrophoresis.
- The DNA is **radiolabeled** to enable **rapid** detection and, in its native state, it migrates **quickly** through the gel matrix.
- Protein binding generates **slower** mobility protein-DNA complexes that resolve as **discrete** bands.


Gel Retardation Assay



Yeast One-hybrid Assay

- One-hybrid screening in yeast is a powerful method to rapidly identify **heterologous transcription factors** that can interact with a **specific** regulatory **DNA sequence** of interest (**bait sequence**).
- In this technique, the interaction between **proteins (prey)** and the **target** DNA sequence (**bait**) is detected by **in vivo reconstitution** of a transcriptional activator that "**turns on**" expression of a **reporter** gene.
- To ensure that DNA binding results in reporter-gene activation, cDNA **expression libraries** are used to produce **hybrids** between the prey and a strong **trans-activating** domain.
- The **advantage** of cloning transcription factors or other DNA-binding proteins by one-hybrid screenings, compared to biochemical techniques, is that the procedure does **not** require specific optimization of **in vitro** conditions.


Yeast One-hybrid Assay




Reporter Gene Assay

- Reporter gene assay – an assay used to study **regulatory promoter** and **enhancer sequences** as well as **transcription factors**.
- Regulatory sequences of interest are **cloned** into a reporter gene-containing vector and are introduced into cells.
- Rationale is that the **level** of reporter gene expression will **correlate** with the **transcriptional activity** of the **DNA fragment of interest**.
- Expression of reporter genes can be measured by:
 1. **Enzyme activity** assay of the expressed enzyme encoded by the reporter gene using chromo-, fluoro- or luminogenic substrates.
 2. **Immunological** assay of the expressed protein encoded by the reporter gene (reporter gene ELISA).
 3. **Histochemical staining** of cells or tissues typically to localize enzymatic activity ectopically expressed from reporter gene constructs in transformed cells.

Reporter Gene Assay



Reporter Genes




- Commonly used reporter genes that induce visually identifiable characteristics usually involve fluorescent and luminescent proteins.
- Examples include:
 - ✓ the gene that encodes jellyfish **green fluorescent protein (GFP)**, which causes cells that express it to glow green under blue light
 - ✓ the enzyme **luciferase**, which catalyzes a reaction with luciferin to produce light
 - ✓ the **red fluorescent protein** from the gene dsRed.

Western Blotting


- Proteins are separated according to **size** by SDS-polyacrylamide gel electrophoresis and transferred from the gel to a filter.
- The filter is incubated with an **antibody** directed against a **protein of interest**.
- The antibody bound to the filter can then be detected by reaction with **various reagents**, such as a radioactive probe that binds to the antibody.

Western Blotting


I. Protein Separation



II. Immunoblotting



III. Detection



South Western Blotting

- South Western blotting is a method for **rapid characterization** of both **DNA-binding proteins** and their **specific sites on genomic DNA**.
- Proteins are separated on a SDS-polyacrylamide gel, renatured by removing SDS in the presence of urea, and blotted onto nitrocellulose by diffusion.
- The genomic DNA region of interest is digested by restriction enzymes selected to produce fragments of **appropriate but different sizes**, which are subsequently **end-labeled** and allowed to bind to the **separated** proteins.
- The specifically-bound DNA is **eluted** from each individual protein-DNA complex and analyzed by acrylamide gel electrophoresis.
- Sequence-specific binding allows the **purification** of the corresponding selectively-bound DNA fragments and **improves** protein-mediated cloning of DNA regulatory sequences.